

Advanced IC

Indwelling and suprapubic vs IC

Patients experiencing urinary retention are faced with a problem. Urine that is left behind in the bladder may contain pathogenic bacteria that can potentially thrive. This may increase the risk of UTIs and can possibly have a negative impact on quality of life for patients.¹⁻³

Intermittent catheterisation (IC) can be used to effectively drain urine by inserting the catheter through the urethra into the bladder. The catheter is inserted only when needed and removed once catheterisation is completed.



“Currently intermittent catheterisation is the recommended treatment for voiding.”

Elena | Senior nurse

There are different bladder management options available for people with urinary retention

1. Suprapubic indwelling catheters

A catheter is inserted into the bladder via the anterior abdominal wall. The catheter is held in place by a balloon and is connected to a urine bag, often placed on the leg.

2. Botox

Often used in combination with other treatment options such as IC. Botox can offer control of the bladder and requires regular repeated injections.

3. Sacral nerve stimulation

A medical treatment that involves the use of electrical impulses to modulate the sacral nerves' function, which controls the bladder, bowel, and pelvic floor muscles.



4. Intermittent catheterisation

Used to effectively drain urine Padsby inserting the catheter through the urethra into the bladder. The catheter is inserted only when needed and removed once the catheterisation is completed.

5. Urethral indwelling catheters

A catheter that's inserted with a balloon through the urethra for bladder drainage. A urine bag will collect the urine.

6. Pads

Often used in combination with other treatment options. This solution may cause leakage, odour, skin issues and a wet feeling, however, it can be an easy solution to manage for users.

7. Medication

Medication will often be used in combination with other treatment options. It can be a gateway to patient acceptance for some users, however for other users, it can have low efficacy.

IC is the gold standard in bladder management⁴

Guidelines

Guidelines from professional associations in Europe and the US promote IC as standard treatment (gold standard) for patients who are unable to empty their bladder.⁴⁻⁶

How

The European Association of Urology Nurses Guidelines state that IC can have a positive impact on patients' quality of life.⁴

Who

The European Association of Urologists (EAU) Guidelines on Neuro-Urology recommend IC for patients who are unable to empty their bladder.⁵

IC has a lower risk of UTIs

Compared with indwelling catheters, IC is associated with a lower risk of UTIs. This has been reported in multiple studies with various patient populations.⁷⁻¹²


Patients using IC show significantly less episodes of bacteriuria and UTIs than indwelling catheters (see Figure 1).¹³ A retrospective study found the rate of symptomatic UTI was lower in patients using IC (70.5%) compared with patients with indwelling catheters (83.3%).¹¹

The same study found that the rate of recurrent UTIs for people using IC was 31.2% whereas the rate in people with an indwelling catheter was 50%.¹¹



Figure 1

Self-catheterisation using IC showed a significantly lower number of bacteriuria and UTIs than indwelling catheters in SCI (spinal cord injury) patients.¹³




IC may reduce the use of antibiotics

Antibiotics are a type of antimicrobial and overuse can lead to Antimicrobial Resistance (AMR). AMR occurs when bacteria, viruses, fungi and parasites change over time and no longer respond to medicines, making infections harder to treat and increasing the risk of disease spread, severe illness and death.¹⁴

Reducing the risk of UTI with IC can help in the global fight against antibiotic resistance. With the high frequency of UTI in patients^{11,15} and antibiotics required for some cases,^{4,16} the consequences can be great.

Antibiotic resistance in bacteria that cause UTI is widespread and a major problem.^{17,18} Reducing the prevalence of UTI can help the antibiotic resistance problem.



IC has a lower risk of vesicoureteral reflux

Vesicoureteral reflux is when urine moves backwards from the bladder to the kidneys. In a retrospective study of spinal cord injury (SCI) patients, the occurrence of vesicoureteral reflux was significantly lower with IC compared with urethral indwelling catheters ($p < 0.001$) and suprapubic indwelling catheters ($p = 0.003$).⁷

IC has a lower risk of bladder stones

Results from a long-term study of SCI patients found that bladder stones were observed in 11% of suprapubic catheter patients, in 6.6% of indwelling catheter patients, and 2% of IC patients.¹⁹

In the same study, a shorter period to bladder stone development was found in patients with an indwelling catheter or suprapubic catheter. The longest time interval was found in patients performing IC.¹⁹

In a retrospective study of patients, bladder stones were significantly less frequent in IC compared with urethral and suprapubic indwelling catheters ($p < 0.001$).⁷

IC can increase patient quality of life

Self-intermittent catheterisation has a better satisfaction level or overall quality of life (QoL) than an indwelling catheter.²⁰

A recent study reported adherence to an intermittent catheterisation (IC) routine of up to 84% after one year and significantly improved QoL in patients performing IC compared with indwelling catheterisation.⁴

IC increases patient independence

One of the most important factors contributing to compliance with long-term IC are continence and the ability to perform IC independently.²¹

References

1. Kennelly, M., et al (2019). Adult Neurogenic Lower Urinary Tract Dysfunction and Intermittent Catheterisation in a Community Setting: Risk Factors Model for Urinary Tract Infections. *Adv Urol*, 2019, 2757862. <https://doi.org/10.1155/2019/2757862>
2. Vasudeva, P, and H Madersbacher. "Factors Implicated in Pathogenesis of Urinary Tract Infections in Neurogenic Bladders: Some Revered, Few Forgotten, Others Ignored." *Neurourol Urodyn*, vol. 33, no. 1, 2014, pp. 95–100.
3. Barber, A. E., et al. "Urinary Tract Infections: Current and Emerging Management Strategies." *Clinical Infectious Diseases*, vol. 57, no. 5, 3 May 2013, pp. 719–724, academic.oup.com/cid/article/57/5/719/310850, <https://doi.org/10.1093/cid/cit284>.
4. Vahr, S. C., et al. Evidence-based Guidelines for Best Practice in Urological Health Care: Catheterisation Urethral Intermittent in Adults, Dilatation Urethral Intermittent in Adults. *EAUN Guideline*, (2024).
5. Skolarikos, A., et al. EAU guidelines on urolithiasis. *EAU Guideline* (2024)
6. Ginsberg, A. D., et al. The AUA/SUFU Guideline on Adult Neurogenic Lower Urinary Tract Dysfunction. *AUA Guideline* (2021)
7. Weld, Kyle J. & Roger R. Dmochowski. "Effect Of Bladder Management On Urological Complications In Spinal Cord Injured Patients." *Journal Of Urology*, Vol. 163, No. 3, Mar. 2000, Pp. 768–772, [https://doi.org/10.1016/S0022-5347\(05\)67800-7](https://doi.org/10.1016/S0022-5347(05)67800-7).
8. Esclar.n De Ruz A, Garc.a Leoni E, Herruzo Cabrera R. Epidemiology and risk factors for urinary tract infection in patients with spinal cord injury. *J Urol* 2000;164(4):1285-89.
9. Wyndaele JJ, De Sy WA, Claessens H. Evaluation of different methods of bladder drainage used in the early care of spinal cord injury patients. *Paraplegia* 1985;23(1):18-26.
10. McGuire EJ, Savastano J. Comparative urological outcome in women with spinal cord injury. *J Urol* 1986;135(4):730-1.
11. Krebs, J et al. "Risk factors for symptomatic urinary tract infections in individuals with chronic neurogenic lower urinary tract dysfunction." *Spinal cord* vol. 54,9 (2016): 682-6. doi:10.1038/sc.2015.214
12. Turi MH, Hanif S, Fasih Q, et al. Proportion of complications in patients practicing clean intermittent self-catheterization (CISC) vs indwelling catheter. *J Pak Med Assoc* 2006;56(9):401-4.
13. Dedeić-Ljubović, Amela, and Mirsada Hukić. "Catheter-related urinary tract infection in patients suffering from spinal cord injuries." *Bosnian journal of basic medical sciences* vol. 9,1 (2009): 2-9. doi:10.17305/bjbm.2009.2849
14. World Health Organization. *Adherence to Long-Term Therapies : Evidence for Action*. Geneva World Health Organization, 2003.
15. Biering-Sørensen, F et al. "Urological situation five years after spinal cord injury." *Scandinavian journal of urology and nephrology* vol. 33,3 (1999): 157-61. doi:10.1080/003655999750015925
16. Pannek J, Blok B, Castro-Diaz D, et al. EAU Guidelines on neurogenic lower urinary tract dysfunction, 2013. Available at: http://uroweb.org/wp-content/uploads/20_Neurogenic-LUTD_LR.pdf [Accessed 13 April 2016].
17. Zowawi, Hosam M et al. "The emerging threat of multidrug-resistant Gram-negative bacteria in urology." *Nature reviews. Urology* vol. 12,10 (2015): 570-84. doi:10.1038/nrurol.2015.199
18. World Health Organisation. "Antimicrobial Resistance." *Who.int*, World Health Organization: WHO, 15 Feb. 2018, www.who.int/en/news-room/fact-sheets/detail/antimicrobial-resistance.
19. Bartel, P., et al. "Bladder stones in patients with spinal cord injury: a long-term study." *Spinal Cord* 52.4 (2014): 295-297.
20. Gharbi, M et al. "Quality of life in neurogenic patients based on different bladder management methods: A review." *Progres en urologie : journal de l'Association francaise d'urologie et de la Societe francaise d'urologie* vol. 32,11 (2022): 784-808. doi:10.1016/j.purol.2022.07.004
21. Perrouin-Verbe, B et al. "Clean intermittent catheterisation from the acute period in spinal cord injury patients. Long term evaluation of urethral and genital tolerance." *Paraplegia* vol. 33,11 (1995): 619-24. doi:10.1038/sc.1995.131

Coloplast A/S, Holtevej 1
3050 Humlebaek, Denmark

Coloplast and the Coloplast logo are trademarks of Coloplast A/S.
© YYYY-MM Coloplast A/S. All rights reserved. PM-37288